



Krabat Runner

- Can be used both indoor and outdoor
- Easy to use for child and caregiver
- Great stability due to three wheels
- Hand brake for child to use
- Combined kick bike and walking bike
- Detachable seat and push bar
- Foldable and height adjustable steering bar
- Height and depth adjustable seat

Follow us on facebook.com/krabaten or Instagram [@krabat_as](https://instagram.com/krabat_as)



Product information

RUNNER SIZE 1

Max. user weight: 35 kg

User age: Approx. 3-8 years

Product weight: 7,8 kg

RUNNER SIZE 2

Max. user weight: 70 kg

User age: From 8 years and up

Product weight: 11,3 kg