

Hepro Trolley

HEPRO
A PART OF **AddLife**



Trolley is a height-adjustable four-wheel trolley that can be rolled over carpet edges and indoor thresholds.

The product is suitable for people who have challenges carrying items between different places indoors and people with impaired hand function.

The robust construction makes it easy for you to move the trolley and use it wherever you want. There is a lot of space on the upper table top, which also has a raised edge that prevents objects from falling off when the trolley is moved. In addition, trolley has a removable metal basket for extra storage space. Lockable wheels for increased stability.

The trolley is easy to assemble, disassemble, clean and adjust to the desired height. The product is 100% recyclable and thus a good choice for the environment.

Product name	Hepro Trolley
Table top, width	50 cm
Table top, depth	34 cm
Overall Width	57 cm
Total length	45 cm
Height table top, min.	52 cm
Height table top, max.	83 cm
Brakes (yes / no)	yes
Swivel wheels (yes / no)	yes
Total weight	7,5 kg
Intended for children (yes / no)	no
Material in table top	steel
Material in frame / legs	steel
Article no.	35600